

FOCUSED ON PRECISION. BUILT FOR

2025

Thank you for your continued trust and partnership.

As we look ahead to 2026, we're excited to take a moment to reflect on the changes and milestones that shaped 2025.

Your projects make a difference and drive us to keep improving every day.

We're grateful to be on this journey with you and look forward to creating many more exciting projects together.

200+

Global Customers

26

Countries

44

States

10+

Trade Shows

Project Spotlight


Completed custom glovebox projects for diverse applications—including battery, nuclear, welding, and R&D.

[See projects →](#)



Our annual maintenance visit includes:

- Leak Tightness (Pressure) inspection to ensure the primary containment barrier integrity
- Functional Verification of All Glovebox Systems (Airflow, Safety, Interlocks, and more)
- Analysis and Record Verification Check
- Vacuum Pump Condition Check and preventive maintenance
- Visual Check for damage
- Seal Maintenance replacement, grease and cleaning
- Software updates if available
- Refresher training on safe glovebox operation and standard best practices
- Detailed inspection and maintenance report provided after service



- 01 On-site Installation & Training
- 02 Preventive Maintenance & Repairs
- 03 Technical Assistance
- 04 OEM Spare Parts



Maintenance Solutions for your glovebox



Long-Term Protection • Reliable Performance

Maintenance Contract

Choose from **Bronze, Silver, or Gold** levels depending on how much coverage and response support you need.

[Explore Maintenance Contract →](#)

Need Service? 870-452-4455 ext 2

Contact US
1-870-452-4455 ext 2
E: Service@inerttechnology.com
Website: www.inert.com

Innovation & Service

Dedicated to keeping customers supplied with critical parts, we've introduced a new ERP system to streamline inventory and enhance service.

[Browse parts →](#)



Trade Show Attendance

Engaged in trade shows across diverse industries, including additive manufacturing, nuclear, battery and more.

[See events →](#)

Custom Project

From C1D1/C1D2 enclosures to multi-bay systems—engineered to spec with CAD reviews, FAT, and on-site installation.

[Having a special project?→](#)



E3i Group

INERT is a member of the E3i group (formerly known as the Gloveboxes Group), alongside Jacomex and EREA. With this new identity, we aim to broaden our reach, drive greater innovation, and deliver more diverse solutions for industries worldwide.

[What is E3i group→](#)

Stay connected to learn how we support industries in driving efficiency and reliability.

[Contact Us](#)

Contact: sales@inerttechnology.com · 1-978-462-4415
1 Industrial Way, Amesbury, MA 01913

